



# PRODUCT CATALOG

Tissue Expanders,  
Breast Implants and Sizers

EFFECTIVE MARCH 2022

**WARNING:**

- Breast implants are not considered lifetime devices. The longer people have them, the greater the chances are that they will develop complications, some of which will require more surgery.
- Breast implants have been associated with the development of a cancer of the immune system called breast implant-associated anaplastic large cell lymphoma (BIA-ALCL). This cancer occurs more commonly in patients with textured breast implants than smooth implants, although rates are not well defined. Some patients have died from BIA-ALCL.
- Patients receiving breast implants have reported a variety of systemic symptoms such as joint pain, muscle aches, confusion, chronic fatigue, autoimmune diseases and others. Individual patient risk for developing these symptoms has not been well established. Some patients report complete resolution of symptoms when the implants are removed without replacement.

The sale and distribution of Mentor Breast Implant Devices are restricted to users and/or user facilities that provide information to patients about the risks and benefits of the device prior to its use in the form and manner specified in approved labeling to be provided by Mentor Worldwide LLC.



# CONTENTS

<b>MemoryGel® Breast Implants</b> .....	4-12
<b>and Resterilizable Smooth Gel Sizers</b>	
Smooth Round Moderate Plus Profile BOOST Breast Implant	
Smooth Round High Profile BOOST Breast Implant	
Smooth Round Moderate Plus Profile Xtra Breast Implant	
Smooth Round Moderate High Profile Xtra Breast Implant	
Smooth Round High Profile Xtra Breast Implant	
Smooth Round Moderate Classic Profile Breast Implant	
Smooth Round Moderate Plus Profile Breast Implant	
Smooth Round High Profile Breast Implant	
Smooth Round Ultra High Profile Breast Implant	
SILTEX® Microtextured Round Mod Plus Xtra Breast Implant	
SILTEX® Microtextured Round High Profile Xtra Breast Implant	
SILTEX® Microtextured Round Moderate Classic Profile Breast Implant	
SILTEX® Microtextured Round Moderate Plus Profile Breast Implant	
SILTEX® Microtextured Round High Profile Breast Implant	
SILTEX® Microtextured Round Ultra High Profile Breast Implant	
<b>MemoryShape® SILTEX® Microtextured Breast Implants</b> .....	13-15
<b>and Resterilizable Gel Sizers</b>	
MM Profile, Medium Height/Moderate Projection Breast Implant	
MM+ Profile, Medium Height/Moderate Plus Projection Breast Implant	
MH Profile, Medium Height/High Projection Breast Implant	
TM+ Profile, Tall Height/Moderate Plus Projection Breast Implant	
LM+ Profile, Low Height/Moderate Plus Projection Breast Implant	
TH Profile, Tall Height/High Projection Breast Implant	
<b>Saline Breast Implants</b> .....	16-17
<b>and Sterile Saline Breast Sizers</b>	
Smooth Round Moderate Profile Breast Implant	
Smooth Round Moderate Plus Profile Breast Implant	
Smooth Round High Profile Breast Implant	
Round Moderate Profile Breast Sizer	
Round Moderate Plus Profile Breast Sizer	
Round High Profile Breast Sizer	
<b>SPECTRUM® Adjustable Saline Breast Implants</b> .....	18
Smooth Round SPECTRUM® Saline Breast Implant	
<b>Tissue Expanders</b> .....	19-22
<b>ARTOURA® Breast Tissue Expanders</b>	
Smooth and SILTEX® High Profile Breast Tissue Expanders	
Smooth and SILTEX® Ultra High Profile Breast Tissue Expanders	
<b>CPX® 4 Breast Tissue Expanders</b>	
Smooth and SILTEX® Low Height Breast Tissue Expanders	
Smooth and SILTEX® Medium Height Breast Tissue Expanders	
Smooth and SILTEX® Tall Height Breast Tissue Expanders	
<b>Smooth Tissue Expanders</b>	
Smooth Round Expander	
Smooth Elliptical Expander	
Smooth Crescent Expander	
Smooth Rectangle Expander	
<b>Accessories</b> .....	23
MENTOR® Volume Sizing System	
Aseptic Transfer System	
Injection Port Detector	
Injection Domes	
Miscellaneous	
<b>MENTOR® Promise Protection Plan</b> .....	24
<b>and MENTOR® Deflation Protection Program</b>	
<b>Important Safety Information</b> .....	25



## MemoryGel BOOST™ Smooth Breast Implants



### Smooth Round Moderate Plus Profile BOOST

Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer†
130 cc	8.8 cm	3.1 cm	<b>SMPB130</b>	RSZMPB130S
160 cc	9.3 cm	3.3 cm	<b>SMPB160</b>	RSZMPB160S
185 cc	9.8 cm	3.3 cm	<b>SMPB185</b>	RSZMPB185S
210 cc	10.3 cm	3.5 cm	<b>SMPB210</b>	RSZMPB210S
235 cc	10.7 cm	3.6 cm	<b>SMPB235</b>	RSZMPB235S
265 cc	11.2 cm	3.6 cm	<b>SMPB265</b>	RSZMPB265S
290 cc	11.4 cm	3.9 cm	<b>SMPB290</b>	RSZMPB290S
320 cc	11.6 cm	4.1 cm	<b>SMPB320</b>	RSZMPB320S
340 cc	12.0 cm	4.1 cm	<b>SMPB340</b>	RSZMPB340S
390 cc	12.7 cm	4.1 cm	<b>SMPB390</b>	RSZMPB390S
430 cc	12.9 cm	4.4 cm	<b>SMPB430</b>	RSZMPB430S
450 cc	13.2 cm	4.5 cm	<b>SMPB450</b>	RSZMPB450S
480 cc	13.4 cm	4.6 cm	<b>SMPB480</b>	RSZMPB480S
530 cc	13.9 cm	4.7 cm	<b>SMPB530</b>	RSZMPB530S
545 cc	14.1 cm	4.8 cm	<b>SMPB545</b>	RSZMPB545S
590 cc	14.3 cm	5.0 cm	<b>SMPB590</b>	RSZMPB590S
630 cc	14.8 cm	5.0 cm	<b>SMPB630</b>	RSZMPB630S
675 cc	15.1 cm	5.0 cm	<b>SMPB675</b>	RSZMPB675S
730 cc	15.6 cm	5.1 cm	<b>SMPB730</b>	RSZMPB730S

### Smooth Round High Profile BOOST



Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer†
200 cc	9.0 cm	4.4 cm	<b>SHPB200</b>	RSZHPB200S
235 cc	9.5 cm	4.7 cm	<b>SHPB235</b>	RSZHPB235S
260 cc	9.8 cm	4.8 cm	<b>SHPB260</b>	RSZHPB260S
285 cc	10.3 cm	4.9 cm	<b>SHPB285</b>	RSZHPB285S
325 cc	10.5 cm	5.3 cm	<b>SHPB325</b>	RSZHPB325S
370 cc	11.0 cm	5.4 cm	<b>SHPB370</b>	RSZHPB370S
405 cc	11.2 cm	5.6 cm	<b>SHPB405</b>	RSZHPB405S
440 cc	11.6 cm	5.8 cm	<b>SHPB440</b>	RSZHPB440S
480 cc	12.2 cm	5.8 cm	<b>SHPB480</b>	RSZHPB480S
520 cc	12.3 cm	6.0 cm	<b>SHPB520</b>	RSZHPB520S
545 cc	12.5 cm	6.1 cm	<b>SHPB545</b>	RSZHPB545S
585 cc	12.8 cm	6.3 cm	<b>SHPB585</b>	RSZHPB585S
635 cc	13.2 cm	6.4 cm	<b>SHPB635</b>	RSZHPB635S
685 cc	13.7 cm	6.4 cm	<b>SHPB685</b>	RSZHPB685S
740 cc	14.1 cm	6.6 cm	<b>SHPB740</b>	RSZHPB740S
775 cc	14.4 cm	6.6 cm	<b>SHPB775</b>	RSZHPB775S

Individual implant dimensions may vary slightly in products of this type.  
Not all units will conform exactly to the dimensions noted above.

\*Shipped sterile.  
Resterilizable up to 10 times.

## MemoryGel® Xtra Smooth Breast Implants and Resterilizable Gel Sizers



### Smooth Round Moderate Plus Profile Xtra

Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer†
115 cc	8.3 cm	3.1 cm	<b>SMPX115</b>	RSZMPX115S
130 cc	8.9 cm	3.1 cm	<b>SMPX130</b>	RSZMPX130S
160 cc	9.5 cm	3.3 cm	<b>SMPX160</b>	RSZMPX160S
190 cc	10.0 cm	3.4 cm	<b>SMPX190</b>	RSZMPX190S
215 cc	10.4 cm	3.6 cm	<b>SMPX215</b>	RSZMPX215S
240 cc	10.8 cm	3.7 cm	<b>SMPX240</b>	RSZMPX240S
270 cc	11.4 cm	3.8 cm	<b>SMPX270</b>	RSZMPX270S
295 cc	11.5 cm	4.0 cm	<b>SMPX295</b>	RSZMPX295S
325 cc	11.9 cm	4.1 cm	<b>SMPX325</b>	RSZMPX325S
350 cc	12.3 cm	4.1 cm	<b>SMPX350</b>	RSZMPX350S
370 cc	12.6 cm	4.1 cm	<b>SMPX370</b>	RSZMPX370S
405 cc	12.7 cm	4.4 cm	<b>SMPX405</b>	RSZMPX405S
440 cc	13.1 cm	4.5 cm	<b>SMPX440</b>	RSZMPX440S
465 cc	13.3 cm	4.7 cm	<b>SMPX465</b>	RSZMPX465S
490 cc	13.6 cm	4.7 cm	<b>SMPX490</b>	RSZMPX490S
525 cc	13.9 cm	4.7 cm	<b>SMPX525</b>	RSZMPX525S
545 cc	14.0 cm	4.9 cm	<b>SMPX545</b>	RSZMPX545S
560 cc	14.3 cm	4.8 cm	<b>SMPX560</b>	RSZMPX560S
605 cc	14.5 cm	5.1 cm	<b>SMPX605</b>	RSZMPX605S
630 cc	14.7 cm	5.1 cm	<b>SMPX630</b>	RSZMPX630S
645 cc	14.9 cm	5.1 cm	<b>SMPX645</b>	RSZMPX645S
685 cc	15.3 cm	5.2 cm	<b>SMPX685</b>	RSZMPX685S
755 cc	15.7 cm	5.4 cm	<b>SMPX755</b>	RSZMPX755S



### Smooth Round Moderate High Profile Xtra

Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer†
130 cc	8.4 cm	3.5 cm	<b>SMHX130</b>	RSZSMHX130S
150 cc	8.8 cm	3.7 cm	<b>SMHX150</b>	RSZSMHX150S
170 cc	9.1 cm	3.8 cm	<b>SMHX170</b>	RSZSMHX170S
185 cc	9.5 cm	3.9 cm	<b>SMHX185</b>	RSZSMHX185S
210 cc	9.8 cm	4.1 cm	<b>SMHX210</b>	RSZSMHX210S
230 cc	10.1 cm	4.2 cm	<b>SMHX230</b>	RSZSMHX230S
255 cc	10.5 cm	4.3 cm	<b>SMHX255</b>	RSZSMHX255S
280 cc	10.8 cm	4.4 cm	<b>SMHX280</b>	RSZSMHX280S
310 cc	11.1 cm	4.6 cm	<b>SMHX310</b>	RSZSMHX310S
330 cc	11.4 cm	4.7 cm	<b>SMHX330</b>	RSZSMHX330S
350 cc	11.6 cm	4.8 cm	<b>SMHX350</b>	RSZSMHX350S
375 cc	11.9 cm	4.9 cm	<b>SMHX375</b>	RSZSMHX375S
400 cc	12.1 cm	4.9 cm	<b>SMHX400</b>	RSZSMHX400S
430 cc	12.5 cm	5.1 cm	<b>SMHX430</b>	RSZSMHX430S
470 cc	12.8 cm	5.2 cm	<b>SMHX470</b>	RSZSMHX470S
500 cc	13.1 cm	5.3 cm	<b>SMHX500</b>	RSZSMHX500S
545 cc	13.4 cm	5.5 cm	<b>SMHX545</b>	RSZSMHX545S
580 cc	13.8 cm	5.5 cm	<b>SMHX580</b>	RSZSMHX580S
630 cc	14.2 cm	5.7 cm	<b>SMHX630</b>	RSZSMHX630S
700 cc	14.6 cm	5.9 cm	<b>SMHX700</b>	RSZSMHX700S
775 cc	15.2 cm	6.0 cm	<b>SMHX775</b>	RSZSMHX775S

Individual implant dimensions may vary slightly in products of this type.  
Not all units will conform exactly to the dimensions noted above.

\*Shipped sterile.  
Resterilizable up to 10 times.

## MemoryGel® Xtra Smooth Breast Implants and Resterilizable Gel Sizers



### Smooth Round High Profile Xtra

Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer†
150 cc	8.3 cm	4.1 cm	<b>SHPX150</b>	RSZSHPX150S
180 cc	8.8 cm	4.4 cm	<b>SHPX180</b>	RSZSHPX180S
200 cc	9.2 cm	4.5 cm	<b>SHPX200</b>	RSZSHPX200S
235 cc	9.5 cm	4.9 cm	<b>SHPX235</b>	RSZSHPX235S
260 cc	9.9 cm	5.0 cm	<b>SHPX260</b>	RSZSHPX260S
285 cc	10.3 cm	5.0 cm	<b>SHPX285</b>	RSZSHPX285S
335 cc	10.7 cm	5.4 cm	<b>SHPX335</b>	RSZSHPX335S
355 cc	11.0 cm	5.5 cm	<b>SHPX355</b>	RSZSHPX355S
380 cc	11.3 cm	5.5 cm	<b>SHPX380</b>	RSZSHPX380S
415 cc	11.5 cm	5.8 cm	<b>SHPX415</b>	RSZSHPX415S
450 cc	11.9 cm	5.9 cm	<b>SHPX450</b>	RSZSHPX450S
465 cc	12.1 cm	5.8 cm	<b>SHPX465</b>	RSZSHPX465S
490 cc	12.5 cm	5.8 cm	<b>SHPX490</b>	RSZSHPX490S
535 cc	12.7 cm	6.0 cm	<b>SHPX535</b>	RSZSHPX535S
560 cc	12.7 cm	6.3 cm	<b>SHPX560</b>	RSZSHPX560S
595 cc	13.1 cm	6.3 cm	<b>SHPX595</b>	RSZSHPX595S
650 cc	13.5 cm	6.5 cm	<b>SHPX650</b>	RSZSHPX650S
700 cc	14.1 cm	6.5 cm	<b>SHPX700</b>	RSZSHPX700S
755 cc	14.4 cm	6.7 cm	<b>SHPX755</b>	RSZSHPX755S
790 cc	14.8 cm	6.7 cm	<b>SHPX790</b>	RSZSHPX790S

Individual implant dimensions may vary slightly in products of this type.  
Not all units will conform exactly to the dimensions noted above.

\*Shipped sterile.  
Resterilizable up to 10 times.

## MemoryGel® Smooth Breast Implants and Resterilizable Gel Sizers



### Smooth Round Moderate Classic Profile

Volume	Diameter	Projection	Catalog #
130 cc	9.6 cm	2.4 cm	<b>3507130MC</b>
150 cc	10.3 cm	2.5 cm	<b>3507150MC</b>
170 cc	10.5 cm	2.6 cm	<b>3507170MC</b>
190 cc	11.0 cm	2.7 cm	<b>3507190MC</b>
215 cc	11.6 cm	2.8 cm	<b>3507215MC</b>
235 cc	11.8 cm	2.9 cm	<b>3507235MC</b>
255 cc	12.2 cm	2.9 cm	<b>3507255MC</b>
275 cc	12.5 cm	3.0 cm	<b>3507275MC</b>
300 cc	12.8 cm	3.1 cm	<b>3507300MC</b>
320 cc	13.0 cm	3.2 cm	<b>3507320MC</b>
340 cc	13.4 cm	3.2 cm	<b>3507340MC</b>
360 cc	13.6 cm	3.3 cm	<b>3507360MC</b>
385 cc	13.9 cm	3.3 cm	<b>3507385MC</b>
405 cc	14.2 cm	3.4 cm	<b>3507405MC</b>
425 cc	14.6 cm	3.4 cm	<b>3507425MC</b>
445 cc	14.8 cm	3.5 cm	<b>3507445MC</b>
470 cc	15.0 cm	3.5 cm	<b>3507470MC</b>
490 cc	15.2 cm	3.7 cm	<b>3507490MC</b>
510 cc	15.4 cm	3.7 cm	<b>3507510MC</b>
555 cc	15.7 cm	3.8 cm	<b>3507555MC</b>
595 cc	16.1 cm	3.8 cm	<b>3507595MC</b>
640 cc	16.5 cm	3.9 cm	<b>3507640MC</b>
680 cc	17.0 cm	3.9 cm	<b>3507680MC</b>
750 cc	17.5 cm	4.0 cm	<b>3507750MC</b>
800 cc	18.0 cm	4.1 cm	<b>3507800MC</b>

#### Resterilizable Gel Sizer\*

RSZ7130MCS
RSZ7150MCS
RSZ7170MCS
RSZ7190MCS
RSZ7215MCS
RSZ7235MCS
RSZ7255MCS
RSZ7275MCS
RSZ7300MCS
RSZ7320MCS
RSZ7340MCS
RSZ7360MCS
RSZ7385MCS
RSZ7405MCS
RSZ7425MCS
RSZ7445MCS
RSZ7470MCS
RSZ7490MCS
RSZ7510MCS
RSZ7555MCS
RSZ7595MCS
RSZ7640MCS
RSZ7680MCS
RSZ7750MCS
RSZ7800MCS

### Smooth Round Moderate Plus Profile



Volume	Diameter	Projection	Catalog #
100 cc	8.2 cm	2.7 cm	<b>3501001BC</b>
125 cc	8.9 cm	2.8 cm	<b>3501251BC</b>
150 cc	9.5 cm	2.9 cm	<b>3501501BC</b>
175 cc	10.0 cm	3.1 cm	<b>3501751BC</b>
200 cc	10.5 cm	3.2 cm	<b>3502001BC</b>
225 cc	10.9 cm	3.3 cm	<b>3502251BC</b>
250 cc	11.3 cm	3.4 cm	<b>3502501BC</b>
275 cc	11.7 cm	3.5 cm	<b>3502751BC</b>
300 cc	12.0 cm	3.6 cm	<b>3503001BC</b>
325 cc	12.3 cm	3.8 cm	<b>3503251BC</b>
350 cc	12.5 cm	3.9 cm	<b>3503501BC</b>
375 cc	12.8 cm	4.0 cm	<b>3503751BC</b>
400 cc	13.1 cm	4.0 cm	<b>3504001BC</b>
425 cc	13.5 cm	4.1 cm	<b>3504251BC</b>
450 cc	13.6 cm	4.2 cm	<b>3504501BC</b>
475 cc	14.0 cm	4.2 cm	<b>3504751BC</b>
500 cc	14.1 cm	4.3 cm	<b>3505001BC</b>
525 cc	14.5 cm	4.4 cm	<b>3505251BC</b>
550 cc	14.6 cm	4.5 cm	<b>3505501BC</b>
575 cc	14.9 cm	4.6 cm	<b>3505751BC</b>
600 cc	15.0 cm	4.6 cm	<b>3506001BC</b>
650 cc	15.4 cm	4.8 cm	<b>3506501BC</b>
700 cc	15.8 cm	4.9 cm	<b>3507001BC</b>
750 cc	16.3 cm	4.9 cm	<b>3507501BC</b>
800 cc	16.5 cm	5.1 cm	<b>3508001BC</b>

#### Resterilizable Gel Sizer\*

RSZ1001S
RSZ1251S
RSZ1501S
RSZ1751S
RSZ2001S
RSZ2251S
RSZ2501S
RSZ2751S
RSZ3001S
RSZ3251S
RSZ3501S
RSZ3751S
RSZ4001S
RSZ4251S
RSZ4501S
RSZ4751S
RSZ5001S
RSZ5251S
RSZ5501S
RSZ5751S
RSZ6001S
RSZ6501S
RSZ7001S
RSZ7501S
RSZ8001S

Individual implant dimensions may vary slightly in products of this type.  
Not all units will conform exactly to the dimensions noted above.

\*Shipped sterile.  
Resterilizable up to 10 times.



## MemoryGel® Smooth Breast Implants and Resterilizable Gel Sizers



### Smooth Round High Profile

Volume	Diameter	Projection	Catalog #
125 cc	8.3 cm	3.5 cm	<b>3501254BC</b>
150 cc	8.8 cm	3.7 cm	<b>3501504BC</b>
175 cc	9.3 cm	3.9 cm	<b>3501754BC</b>
200 cc	9.7 cm	4.0 cm	<b>3502004BC</b>
225 cc	10.1 cm	4.2 cm	<b>3502254BC</b>
250 cc	10.5 cm	4.3 cm	<b>3502504BC</b>
275 cc	10.8 cm	4.4 cm	<b>3502754BC</b>
300 cc	11.1 cm	4.5 cm	<b>3503004BC</b>
325 cc	11.4 cm	4.6 cm	<b>3503254BC</b>
350 cc	11.7 cm	4.8 cm	<b>3503504BC</b>
375 cc	12.0 cm	4.8 cm	<b>3503754BC</b>
400 cc	12.2 cm	5.0 cm	<b>3504004BC</b>
425 cc	12.5 cm	5.0 cm	<b>3504254BC</b>
450 cc	12.8 cm	5.1 cm	<b>3504504BC</b>
475 cc	12.9 cm	5.3 cm	<b>3504754BC</b>
500 cc	13.2 cm	5.3 cm	<b>3505004BC</b>
550 cc	13.6 cm	5.5 cm	<b>3505504BC</b>
600 cc	14.0 cm	5.6 cm	<b>3506004BC</b>
650 cc	14.4 cm	5.7 cm	<b>3506504BC</b>
700 cc	14.8 cm	5.8 cm	<b>3507004BC</b>
750 cc	15.1 cm	5.9 cm	<b>3507504BC</b>
800 cc	15.5 cm	6.0 cm	<b>3508004BC</b>

#### Resterilizable Gel Sizer\*

RSZ1254S
RSZ1504S
RSZ1754S
RSZ2004S
RSZ2254S
RSZ2504S
RSZ2754S
RSZ3004S
RSZ3254S
RSZ3504S
RSZ3754S
RSZ4004S
RSZ4254S
RSZ4504S
RSZ4754S
RSZ5004S
RSZ5504S
RSZ6004S
RSZ6504S
RSZ7004S
RSZ7504S
RSZ8004S

### Smooth Round Ultra High Profile



Volume	Diameter	Projection	Catalog #
135 cc	7.8 cm	4.1 cm	<b>3505135BC</b>
160 cc	8.2 cm	4.3 cm	<b>3505160BC</b>
185 cc	8.4 cm	4.4 cm	<b>3505185BC</b>
215 cc	8.7 cm	4.5 cm	<b>3505215BC</b>
240 cc	9.0 cm	4.7 cm	<b>3505240BC</b>
270 cc	9.3 cm	4.8 cm	<b>3505270BC</b>
295 cc	9.7 cm	5.0 cm	<b>3505295BC</b>
320 cc	9.8 cm	5.1 cm	<b>3505320BC</b>
350 cc	10.1 cm	5.2 cm	<b>3505350BC</b>
375 cc	10.4 cm	5.3 cm	<b>3505375BC</b>
400 cc	10.6 cm	5.4 cm	<b>3505400BC</b>
430 cc	10.9 cm	5.6 cm	<b>3505430BC</b>
455 cc	11.2 cm	5.7 cm	<b>3505455BC</b>
480 cc	11.4 cm	5.8 cm	<b>3505480BC</b>
535 cc	12.0 cm	6.1 cm	<b>3505535BC</b>
590 cc	12.5 cm	6.3 cm	<b>3505590BC</b>
650 cc	12.9 cm	6.4 cm	<b>3505650BC</b>
700 cc	13.3 cm	6.5 cm	<b>3505700BC</b>
750 cc	13.7 cm	6.6 cm	<b>3505750BC</b>
800 cc	14.1 cm	6.7 cm	<b>3505800BC</b>

#### Resterilizable Gel Sizer\*

RSZ5135S
RSZ5160S
RSZ5185S
RSZ5215S
RSZ5240S
RSZ5270S
RSZ5295S
RSZ5320S
RSZ5350S
RSZ5375S
RSZ5400S
RSZ5430S
RSZ5455S
RSZ5480S
RSZ5535S
RSZ5590S
RSZ5650S
RSZ5700S
RSZ5750S
RSZ5800S

Individual implant dimensions may vary slightly in products of this type.  
Not all units will conform exactly to the dimensions noted above.

\*Shipped sterile.  
Resterilizable up to 10 times.

## MemoryGel® SILTEX® Microtextured Xtra Breast Implants and Resterilizable Gel Sizers

### SILTEX® Microtextured Round Moderate Plus Xtra



Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer†
115 cc	8.4 cm	3.2 cm	TMPX115	RSZMPX115S
130 cc	9.1 cm	3.1 cm	TMPX130	RSZMPX130S
160 cc	9.5 cm	3.3 cm	TMPX160	RSZMPX160S
190 cc	10.0 cm	3.5 cm	TMPX190	RSZMPX190S
215 cc	10.5 cm	3.6 cm	TMPX215	RSZMPX215S
240 cc	10.9 cm	3.8 cm	TMPX240	RSZMPX240S
270 cc	11.4 cm	3.9 cm	TMPX270	RSZMPX270S
295 cc	11.6 cm	4.1 cm	TMPX295	RSZMPX295S
325 cc	11.9 cm	4.2 cm	TMPX325	RSZMPX325S
350 cc	12.2 cm	4.3 cm	TMPX350	RSZMPX350S
370 cc	12.7 cm	4.3 cm	TMPX370	RSZMPX370S
405 cc	12.8 cm	4.8 cm	TMPX405	RSZMPX405S
440 cc	13.1 cm	4.8 cm	TMPX440	RSZMPX440S
465 cc	13.4 cm	4.8 cm	TMPX465	RSZMPX465S
490 cc	13.6 cm	4.8 cm	TMPX490	RSZMPX490S
525 cc	13.7 cm	5.1 cm	TMPX525	RSZMPX525S
545 cc	14.1 cm	5.0 cm	TMPX545	RSZMPX545S
560 cc	14.3 cm	5.0 cm	TMPX560	RSZMPX560S
605 cc	14.4 cm	5.3 cm	TMPX605	RSZMPX605S
630 cc	14.6 cm	5.4 cm	TMPX630	RSZMPX630S
645 cc	15.0 cm	5.3 cm	TMPX645	RSZMPX645S
685 cc	15.4 cm	5.3 cm	TMPX685	RSZMPX685S
755 cc	15.7 cm	5.6 cm	TMPX755	RSZMPX755S

### SILTEX® Microtextured Round High Profile Xtra



Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer†
150 cc	8.4 cm	4.2 cm	THPX150	RSZTHPX150S
175 cc	9.0 cm	4.4 cm	THPX175	RSZTHPX175S
200 cc	9.4 cm	4.5 cm	THPX200	RSZTHPX200S
230 cc	9.8 cm	4.8 cm	THPX230	RSZTHPX230S
255 cc	10.1 cm	4.9 cm	THPX255	RSZTHPX255S
285 cc	10.5 cm	5.1 cm	THPX285	RSZTHPX285S
325 cc	10.8 cm	5.4 cm	THPX325	RSZTHPX325S
340 cc	11.1 cm	5.4 cm	THPX340	RSZTHPX340S
365 cc	11.4 cm	5.5 cm	THPX365	RSZTHPX365S
405 cc	11.6 cm	5.9 cm	THPX405	RSZTHPX405S
425 cc	11.9 cm	5.8 cm	THPX425	RSZTHPX425S
455 cc	12.1 cm	6.0 cm	THPX455	RSZTHPX455S
470 cc	12.3 cm	5.9 cm	THPX470	RSZTHPX470S
515 cc	12.6 cm	6.2 cm	THPX515	RSZTHPX515S
545 cc	12.8 cm	6.3 cm	THPX545	RSZTHPX545S
570 cc	13.2 cm	6.3 cm	THPX570	RSZTHPX570S
620 cc	13.4 cm	6.5 cm	THPX620	RSZTHPX620S
680 cc	14.0 cm	6.6 cm	THPX680	RSZTHPX680S
725 cc	14.5 cm	6.7 cm	THPX725	RSZTHPX725S
765 cc	14.7 cm	6.7 cm	THPX765	RSZTHPX765S

Individual implant dimensions may vary slightly in products of this type.  
Not all units will conform exactly to the dimensions noted above.

\*Shipped sterile.  
Resterilizable up to 10 times.

## MemoryGel® SILTEX® Microtextured Breast Implants and Resterilizable Gel Sizers

### SILTEX® Microtextured Round Moderate Classic Profile



Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer*
130 cc	9.6 cm	2.5 cm	<b>3541307MC</b>	RSZ7130MCS
150 cc	10.2 cm	2.6 cm	<b>3541507MC</b>	RSZ7150MCS
170 cc	10.6 cm	2.7 cm	<b>3541707MC</b>	RSZ7170MCS
190 cc	11.0 cm	2.8 cm	<b>3541907MC</b>	RSZ7190MCS
215 cc	11.5 cm	2.9 cm	<b>3542157MC</b>	RSZ7215MCS
235 cc	11.7 cm	3.0 cm	<b>3542357MC</b>	RSZ7235MCS
255 cc	12.1 cm	3.1 cm	<b>3542557MC</b>	RSZ7255MCS
275 cc	12.3 cm	3.2 cm	<b>3542757MC</b>	RSZ7275MCS
300 cc	12.6 cm	3.3 cm	<b>3543007MC</b>	RSZ7300MCS
320 cc	12.9 cm	3.3 cm	<b>3543207MC</b>	RSZ7320MCS
340 cc	13.2 cm	3.4 cm	<b>3543407MC</b>	RSZ7340MCS
360 cc	13.4 cm	3.4 cm	<b>3543607MC</b>	RSZ7360MCS
385 cc	13.7 cm	3.5 cm	<b>3543857MC</b>	RSZ7385MCS
405 cc	14.0 cm	3.6 cm	<b>3544057MC</b>	RSZ7405MCS
425 cc	14.3 cm	3.6 cm	<b>3544257MC</b>	RSZ7425MCS
445 cc	14.5 cm	3.7 cm	<b>3544457MC</b>	RSZ7445MCS
470 cc	14.7 cm	3.7 cm	<b>3544707MC</b>	RSZ7470MCS
490 cc	14.9 cm	3.8 cm	<b>3544907MC</b>	RSZ7490MCS
510 cc	15.1 cm	3.8 cm	<b>3545107MC</b>	RSZ7510MCS
555 cc	15.4 cm	4.0 cm	<b>3545557MC</b>	RSZ7555MCS
595 cc	15.8 cm	4.1 cm	<b>3545957MC</b>	RSZ7595MCS
640 cc	16.3 cm	4.2 cm	<b>3546407MC</b>	RSZ7640MCS
680 cc	16.7 cm	4.2 cm	<b>3546807MC</b>	RSZ7680MCS
750 cc	17.2 cm	4.3 cm	<b>3547507MC</b>	RSZ7750MCS
800 cc	17.6 cm	4.4 cm	<b>3548007MC</b>	RSZ7800MCS

### SILTEX® Microtextured Round Moderate Plus Profile



Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer*
100 cc	8.1 cm	2.7 cm	<b>3541001</b>	RSZ1001S
125 cc	8.8 cm	2.9 cm	<b>3541251</b>	RSZ1251S
150 cc	9.4 cm	3.0 cm	<b>3541501</b>	RSZ1501S
175 cc	10.0 cm	3.2 cm	<b>3541751</b>	RSZ1751S
200 cc	10.5 cm	3.3 cm	<b>3542001</b>	RSZ2001S
225 cc	10.9 cm	3.5 cm	<b>3542251</b>	RSZ2251S
250 cc	11.3 cm	3.6 cm	<b>3542501</b>	RSZ2501S
275 cc	11.7 cm	3.7 cm	<b>3542751</b>	RSZ2751S
300 cc	12.0 cm	3.7 cm	<b>3543001</b>	RSZ3001S
325 cc	12.3 cm	3.8 cm	<b>3543251</b>	RSZ3251S
350 cc	12.6 cm	3.8 cm	<b>3543501</b>	RSZ3501S
375 cc	12.9 cm	3.9 cm	<b>3543751</b>	RSZ3751S
400 cc	13.2 cm	4.0 cm	<b>3544001</b>	RSZ4001S
425 cc	13.5 cm	4.2 cm	<b>3544251</b>	RSZ4251S
450 cc	13.7 cm	4.1 cm	<b>3544501</b>	RSZ4501S
475 cc	13.9 cm	4.4 cm	<b>3544751</b>	RSZ4751S
500 cc	14.1 cm	4.2 cm	<b>3545001</b>	RSZ5001S
525 cc	14.4 cm	4.6 cm	<b>3545251</b>	RSZ5251S
550 cc	14.4 cm	4.4 cm	<b>3545501</b>	RSZ5501S
575 cc	14.9 cm	4.7 cm	<b>3545751</b>	RSZ5751S
600 cc	14.7 cm	4.5 cm	<b>3546001</b>	RSZ6001S
650 cc	15.5 cm	4.9 cm	<b>3546501</b>	RSZ6501S
700 cc	15.7 cm	4.8 cm	<b>3547001</b>	RSZ7001S
750 cc	16.2 cm	5.1 cm	<b>3547501</b>	RSZ7501S
800 cc	16.6 cm	5.0 cm	<b>3548001</b>	RSZ8001S

Individual implant dimensions may vary slightly in products of this type.  
Not all units will conform exactly to the dimensions noted above.

\*Shipped sterile.  
Resterilizable up to 10 times.

## MemoryGel® SILTEX® Microtextured Breast Implants and Resterilizable Gel Sizers



### SILTEX® Microtextured Round High Profile

Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer*
125 cc	8.4 cm	3.6 cm	<b>3544125</b>	RSZ1254S
150 cc	8.9 cm	3.8 cm	<b>3544150</b>	RSZ1504S
175 cc	9.4 cm	4.0 cm	<b>3544175</b>	RSZ1754S
200 cc	9.9 cm	4.1 cm	<b>3544200</b>	RSZ2004S
225 cc	10.2 cm	4.3 cm	<b>3544225</b>	RSZ2254S
250 cc	10.5 cm	4.5 cm	<b>3544250</b>	RSZ2504S
275 cc	10.9 cm	4.6 cm	<b>3544275</b>	RSZ2754S
300 cc	11.1 cm	4.7 cm	<b>3544300</b>	RSZ3004S
325 cc	11.5 cm	4.8 cm	<b>3544325</b>	RSZ3254S
350 cc	11.7 cm	4.9 cm	<b>3544350</b>	RSZ3504S
375 cc	12.0 cm	5.0 cm	<b>3544375</b>	RSZ3754S
400 cc	12.3 cm	5.1 cm	<b>3544400</b>	RSZ4004S
425 cc	12.5 cm	5.2 cm	<b>3544425</b>	RSZ4254S
450 cc	12.7 cm	5.2 cm	<b>3544450</b>	RSZ4504S
475 cc	12.9 cm	5.4 cm	<b>3544475</b>	RSZ4754S
500 cc	13.2 cm	5.4 cm	<b>3544500</b>	RSZ5004S
550 cc	13.5 cm	5.6 cm	<b>3544550</b>	RSZ5504S
600 cc	14.0 cm	5.7 cm	<b>3544600</b>	RSZ6004S
650 cc	14.3 cm	5.8 cm	<b>3544650</b>	RSZ6504S
700 cc	14.7 cm	6.0 cm	<b>3544700</b>	RSZ7004S
750 cc	15.1 cm	6.0 cm	<b>3544750</b>	RSZ7504S
800 cc	15.4 cm	6.3 cm	<b>3544800</b>	RSZ8004S



### SILTEX® Microtextured Round Ultra High Profile

Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer*
135 cc	8.0 cm	4.3 cm	<b>3545135</b>	RSZ5135S
160 cc	8.4 cm	4.5 cm	<b>3545160</b>	RSZ5160S
185 cc	8.6 cm	4.6 cm	<b>3545185</b>	RSZ5185S
215 cc	8.9 cm	4.7 cm	<b>3545215</b>	RSZ5215S
240 cc	9.2 cm	4.9 cm	<b>3545240</b>	RSZ5240S
270 cc	9.5 cm	5.0 cm	<b>3545270</b>	RSZ5270S
295 cc	9.9 cm	5.2 cm	<b>3545295</b>	RSZ5295S
320 cc	10.0 cm	5.3 cm	<b>3545320</b>	RSZ5320S
350 cc	10.3 cm	5.4 cm	<b>3545350</b>	RSZ5350S
375 cc	10.5 cm	5.5 cm	<b>3545375</b>	RSZ5375S
400 cc	10.8 cm	5.6 cm	<b>3545400</b>	RSZ5400S
430 cc	11.1 cm	5.8 cm	<b>3545430</b>	RSZ5430S
455 cc	11.3 cm	5.9 cm	<b>3545455</b>	RSZ5455S
480 cc	11.6 cm	6.0 cm	<b>3545480</b>	RSZ5480S
535 cc	12.1 cm	6.3 cm	<b>3545535</b>	RSZ5535S
590 cc	12.6 cm	6.5 cm	<b>3545590</b>	RSZ5590S
650 cc	13.0 cm	6.6 cm	<b>3545650</b>	RSZ5650S
700 cc	13.4 cm	6.7 cm	<b>3545700</b>	RSZ5700S

Individual implant dimensions may vary slightly in products of this type.  
Not all units will conform exactly to the dimensions noted above.

\*Shipped sterile.  
Resterilizable up to 10 times.

## MemoryShape® SILTEX® Microtextured Breast Implants and Resterilizable Gel Sizers



### MM Profile, Medium Height/Moderate Projection

Vol.	Width	Height	Proj.	Catalog #	Resterilizable Gel Sizer*
120 cc	9.0 cm	8.5 cm	3.3 cm	<b>3540908</b>	RSZ0908S
135 cc	9.5 cm	8.9 cm	3.5 cm	<b>3540958</b>	RSZ0958S
155 cc	10.0 cm	9.4 cm	3.7 cm	<b>3541008</b>	RSZ1008S
180 cc	10.5 cm	9.9 cm	3.8 cm	<b>3541058</b>	RSZ1058S
215 cc	11.0 cm	10.3 cm	3.9 cm	<b>3541108</b>	RSZ1108S
245 cc	11.5 cm	10.8 cm	4.0 cm	<b>3541158</b>	RSZ1158S
280 cc	12.0 cm	11.3 cm	4.2 cm	<b>3541208</b>	RSZ1208S
315 cc	12.5 cm	11.8 cm	4.4 cm	<b>3541258</b>	RSZ1258S
355 cc	13.0 cm	12.2 cm	4.6 cm	<b>3541308</b>	RSZ1308S
395 cc	13.5 cm	12.7 cm	4.7 cm	<b>3541358</b>	RSZ1358S
440 cc	14.0 cm	13.2 cm	4.9 cm	<b>3541408</b>	RSZ1408S
480 cc	14.5 cm	13.6 cm	5.0 cm	<b>3541458</b>	RSZ1458S
530 cc	15.0 cm	14.1 cm	5.2 cm	<b>3541508</b>	RSZ1508S
640 cc	16.0 cm	15.0 cm	5.6 cm	<b>3541608</b>	RSZ1608S
775 cc	17.0 cm	16.0 cm	5.9 cm	<b>3541708</b>	RSZ1708S

### MM+ Profile, Medium Height/ Moderate Plus Projection



Vol.	Width	Height	Proj.	Catalog #	Resterilizable Gel Sizer*
140 cc	9.0 cm	8.5 cm	3.8 cm	<b>3340905</b>	RSZ0905S
165 cc	9.5 cm	8.9 cm	4.0 cm	<b>3340955</b>	RSZ0955S
195 cc	10.0 cm	9.4 cm	4.2 cm	<b>3341005</b>	RSZ1005S
225 cc	10.5 cm	9.9 cm	4.4 cm	<b>3341055</b>	RSZ1055S
255 cc	11.0 cm	10.3 cm	4.7 cm	<b>3341105</b>	RSZ1105S
295 cc	11.5 cm	10.8 cm	4.9 cm	<b>3341155</b>	RSZ1155S
330 cc	12.0 cm	11.3 cm	5.1 cm	<b>3341205</b>	RSZ1205S
375 cc	12.5 cm	11.8 cm	5.3 cm	<b>3341255</b>	RSZ1255S
420 cc	13.0 cm	12.2 cm	5.5 cm	<b>3341305</b>	RSZ1305S
475 cc	13.5 cm	12.7 cm	5.7 cm	<b>3341355</b>	RSZ1355S
525 cc	14.0 cm	13.2 cm	5.9 cm	<b>3341405</b>	RSZ1405S
585 cc	14.5 cm	13.6 cm	6.1 cm	<b>3341455</b>	RSZ1455S
650 cc	15.0 cm	14.1 cm	6.3 cm	<b>3341505</b>	RSZ1505S

Individual implant dimensions may vary slightly in products of this type.  
Not all units will conform exactly to the dimensions noted above.

\*Shipped sterile.  
Resterilizable up to 10 times.

## MemoryShape® SILTEX® Microtextured Breast Implants and Resterilizable Gel Sizers



### MH Profile, Medium Height/High Projection

Vol.	Width	Height	Proj.	Catalog #	Resterilizable Gel Sizer*
165 cc	9.0 cm	8.5 cm	4.6 cm	<b>3340902</b>	RSZ0902S
195 cc	9.5 cm	8.9 cm	4.8 cm	<b>3340952</b>	RSZ0952S
225 cc	10.0 cm	9.4 cm	5.1 cm	<b>3341002</b>	RSZ1002S
260 cc	10.5 cm	9.9 cm	5.3 cm	<b>3341052</b>	RSZ1052S
300 cc	11.0 cm	10.3 cm	5.6 cm	<b>3341102</b>	RSZ1102S
345 cc	11.5 cm	10.8 cm	5.8 cm	<b>3341152</b>	RSZ1152S
390 cc	12.0 cm	11.3 cm	6.0 cm	<b>3341202</b>	RSZ1202S
440 cc	12.5 cm	11.8 cm	6.2 cm	<b>3341252</b>	RSZ1252S
495 cc	13.0 cm	12.2 cm	6.5 cm	<b>3341302</b>	RSZ1302S
555 cc	13.5 cm	12.7 cm	6.7 cm	<b>3341352</b>	RSZ1352S
620 cc	14.0 cm	13.2 cm	6.9 cm	<b>3341402</b>	RSZ1402S
685 cc	14.5 cm	13.6 cm	7.1 cm	<b>3341452</b>	RSZ1452S



### TH Profile, Tall Height/High Projection

Vol.	Width	Height	Proj.	Catalog #	Resterilizable Gel Sizer*
180 cc	9.0 cm	9.4 cm	4.6 cm	<b>3340904</b>	RSZ4090S
215 cc	9.5 cm	9.9 cm	4.8 cm	<b>3340954</b>	RSZ4095S
250 cc	10.0 cm	10.4 cm	5.1 cm	<b>3341004</b>	RSZ4100S
290 cc	10.5 cm	10.9 cm	5.3 cm	<b>3341054</b>	RSZ4105S
330 cc	11.0 cm	11.5 cm	5.6 cm	<b>3341104</b>	RSZ4110S
380 cc	11.5 cm	12.0 cm	5.8 cm	<b>3341154</b>	RSZ4115S
430 cc	12.0 cm	12.5 cm	6.0 cm	<b>3341204</b>	RSZ4120S
485 cc	12.5 cm	13.0 cm	6.2 cm	<b>3341254</b>	RSZ4125S
545 cc	13.0 cm	13.5 cm	6.5 cm	<b>3341304</b>	RSZ4130S
610 cc	13.5 cm	14.1 cm	6.7 cm	<b>3341354</b>	RSZ4135S
680 cc	14.0 cm	14.6 cm	6.9 cm	<b>3341404</b>	RSZ4140S
755 cc	14.5 cm	15.1 cm	7.1 cm	<b>3341454</b>	RSZ4145S

Individual implant dimensions may vary slightly in products of this type.  
Not all units will conform exactly to the dimensions noted above.

\*Shipped sterile.  
Resterilizable up to 10 times.

## MemoryShape® SILTEX® Microtextured Breast Implants and Resterilizable Gel Sizers



### TM+ Profile, Tall Height/Moderate Plus Projection

Vol.	Width	Height	Proj.	Catalog #	Resterilizable Gel Sizer*
145 cc	9.0 cm	9.4 cm	3.8 cm	<b>3340909</b>	RSZ0909S
175 cc	9.5 cm	9.9 cm	4.0 cm	<b>3340959</b>	RSZ0959S
205 cc	10.0 cm	10.4 cm	4.2 cm	<b>3341009</b>	RSZ1009S
235 cc	10.5 cm	10.9 cm	4.4 cm	<b>3341059</b>	RSZ1059S
270 cc	11.0 cm	11.5 cm	4.7 cm	<b>3341109</b>	RSZ1109S
305 cc	11.5 cm	12.0 cm	4.9 cm	<b>3341159</b>	RSZ1159S
350 cc	12.0 cm	12.5 cm	5.1 cm	<b>3341209</b>	RSZ1209S
395 cc	12.5 cm	13.0 cm	5.3 cm	<b>3341259</b>	RSZ1259S
445 cc	13.0 cm	13.5 cm	5.5 cm	<b>3341309</b>	RSZ1309S
495 cc	13.5 cm	14.1 cm	5.7 cm	<b>3341359</b>	RSZ1359S
555 cc	14.0 cm	14.6 cm	5.9 cm	<b>3341409</b>	RSZ1409S
615 cc	14.5 cm	15.1 cm	6.1 cm	<b>3341459</b>	RSZ1459S
680 cc	15.0 cm	15.6 cm	6.3 cm	<b>3341509</b>	RSZ1509S

### LM+ Profile, Low Height/Moderate Plus Projection



Vol.	Width	Height	Proj.	Catalog #	Resterilizable Gel Sizer*
125 cc	9.0 cm	8.0 cm	3.8 cm	<b>3340907</b>	RSZ0907S
145 cc	9.5 cm	8.4 cm	4.0 cm	<b>3340957</b>	RSZ0957S
170 cc	10.0 cm	8.8 cm	4.2 cm	<b>3341007</b>	RSZ1007S
195 cc	10.5 cm	9.3 cm	4.4 cm	<b>3341057</b>	RSZ1057S
225 cc	11.0 cm	9.7 cm	4.7 cm	<b>3341107</b>	RSZ1107S
255 cc	11.5 cm	10.2 cm	4.9 cm	<b>3341157</b>	RSZ1157S
290 cc	12.0 cm	10.6 cm	5.1 cm	<b>3341207</b>	RSZ1207S
330 cc	12.5 cm	11.1 cm	5.3 cm	<b>3341257</b>	RSZ1257S
370 cc	13.0 cm	11.5 cm	5.5 cm	<b>3341307</b>	RSZ1307S
415 cc	13.5 cm	12.0 cm	5.7 cm	<b>3341357</b>	RSZ1357S
465 cc	14.0 cm	12.4 cm	5.9 cm	<b>3341407</b>	RSZ1407S
515 cc	14.5 cm	12.8 cm	6.1 cm	<b>3341457</b>	RSZ1457S
570 cc	15.0 cm	13.3 cm	6.3 cm	<b>3341507</b>	RSZ1507S
690 cc	16.0 cm	14.2 cm	6.8 cm	<b>3341607</b>	RSZ1607S

Individual implant dimensions may vary slightly in products of this type.  
Not all units will conform exactly to the dimensions noted above.

\*Shipped sterile.  
Resterilizable up to 10 times.

## SALINE BREAST IMPLANTS

### Smooth Round Moderate Profile



NOMINAL FILL			MAXIMUM FILL			Catalog #	Sterile Saline Sizer†
Volume	Diameter	Proj.	Vol.	Diam.	Proj.		
125 cc	9.5 cm	3.0 cm	150 cc	9.1 cm	3.4 cm	<b>3501610</b>	351-125SZ
150 cc	10.0 cm	3.1 cm	175 cc	9.6 cm	3.7 cm	<b>3501615</b>	
175 cc	10.6 cm	3.3 cm	200 cc	10.0 cm	3.9 cm	<b>3501620</b>	351-175SZ
200 cc	11.0 cm	3.4 cm	225 cc	10.4 cm	3.9 cm	<b>3501625</b>	
225 cc	11.5 cm	3.5 cm	250 cc	10.9 cm	4.0 cm	<b>3501630</b>	351-225SZ
250 cc	11.9 cm	3.6 cm	275 cc	11.2 cm	4.1 cm	<b>3501635</b>	
275 cc	12.3 cm	3.7 cm	300 cc	11.7 cm	4.2 cm	<b>3501640</b>	351-275SZ
300 cc	12.6 cm	3.7 cm	325 cc	12.1 cm	4.3 cm	<b>3501645</b>	
325 cc	13.0 cm	3.8 cm	375 cc	12.3 cm	4.5 cm	<b>3501650</b>	351-325SZ
350 cc	13.3 cm	3.9 cm	400 cc	12.6 cm	4.7 cm	<b>3501655</b>	
375 cc	13.6 cm	4.0 cm	425 cc	12.8 cm	4.8 cm	<b>3501660</b>	351-375SZ
425 cc	14.2 cm	4.1 cm	475 cc	13.5 cm	4.9 cm	<b>3501670</b>	351-425SZ
475 cc	14.8 cm	4.2 cm	525 cc	14.2 cm	4.9 cm	<b>3501680</b>	351-475SZ
525 cc	15.0 cm	4.2 cm	575 cc	14.5 cm	5.1 cm	<b>3501685</b>	
575 cc	15.0 cm	4.5 cm	625 cc	14.8 cm	5.1 cm	<b>3501690</b>	
625 cc	15.2 cm	4.6 cm	700 cc	14.9 cm	5.6 cm	<b>3501695</b>	
700 cc	15.6 cm	4.9 cm	775 cc	14.9 cm	5.8 cm	<b>3501697</b>	

†Single Use Only

### Smooth Round Moderate Plus Profile



NOMINAL FILL			MAXIMUM FILL			Catalog #	Sterile Saline Sizer†
Volume	Diameter	Proj.	Vol.	Diam.	Proj.		
175 cc	9.5 cm	3.5 cm	210 cc	9.3 cm	4.2 cm	<b>3502175</b>	351-2175SZ
200 cc	10.0 cm	3.7 cm	240 cc	9.7 cm	4.4 cm	<b>3502200</b>	351-2200SZ
225 cc	10.4 cm	3.8 cm	270 cc	10.1 cm	4.6 cm	<b>3502225</b>	351-2225SZ
250 cc	10.8 cm	4.0 cm	300 cc	10.5 cm	4.7 cm	<b>3502250</b>	351-2250SZ
275 cc	11.0 cm	4.1 cm	330 cc	10.8 cm	4.9 cm	<b>3502275</b>	351-2275SZ
300 cc	11.5 cm	4.3 cm	360 cc	11.2 cm	5.0 cm	<b>3502300</b>	351-2300SZ
325 cc	11.9 cm	4.4 cm	390 cc	11.5 cm	5.2 cm	<b>3502325</b>	351-2325SZ
350 cc	12.1 cm	4.5 cm	420 cc	11.7 cm	5.3 cm	<b>3502350</b>	351-2350SZ
375 cc	12.3 cm	4.6 cm	450 cc	12.0 cm	5.4 cm	<b>3502375</b>	351-2375SZ
400 cc	12.6 cm	4.7 cm	480 cc	12.3 cm	5.5 cm	<b>3502400</b>	351-2400SZ
425 cc	12.9 cm	4.8 cm	510 cc	12.5 cm	5.6 cm	<b>3502425</b>	351-2425SZ
450 cc	13.0 cm	4.9 cm	540 cc	12.8 cm	5.7 cm	<b>3502450</b>	351-2450SZ
475 cc	13.3 cm	5.0 cm	570 cc	13.0 cm	5.8 cm	<b>3502475</b>	351-2475SZ
500 cc	13.6 cm	5.1 cm	600 cc	13.2 cm	5.9 cm	<b>3502500</b>	351-2500SZ
550 cc	14.0 cm	5.3 cm	660 cc	13.7 cm	6.1 cm	<b>3502550</b>	351-2550SZ
600 cc	14.5 cm	5.5 cm	720 cc	14.1 cm	6.3 cm	<b>3502600</b>	351-2600SZ
650 cc	14.8 cm	5.6 cm	780 cc	14.4 cm	6.5 cm	<b>3502650</b>	351-2650SZ
700 cc	15.2 cm	5.8 cm	840 cc	14.8 cm	6.6 cm	<b>3502700</b>	351-2700SZ
750 cc	15.6 cm	5.9 cm	900 cc	15.1 cm	6.8 cm	<b>3502750</b>	351-2750SZ
800 cc	16.1 cm	6.0 cm	960 cc	15.5 cm	7.0 cm	<b>3502800</b>	351-2800SZ

Individual implant dimensions may vary slightly in products of this type. Not all units will conform exactly to the dimensions noted above.

†Single Use Only



## SALINE BREAST IMPLANTS

### Smooth Round High Profile




NOMINAL FILL			MAXIMUM FILL			Catalog #	Sterile Saline Sizer†
Volume	Diameter	Proj.	Vol.	Diam.	Proj.		
170 cc	9.0 cm	4.0 cm	200 cc	8.8 cm	4.9 cm	<b>3503170</b>	
190 cc	9.3 cm	4.1 cm	225 cc	9.1 cm	5.1 cm	<b>3503190</b>	351-3190SZ
210 cc	9.6 cm	4.2 cm	250 cc	9.5 cm	5.3 cm	<b>3503210</b>	
230 cc	10.0 cm	4.3 cm	275 cc	9.8 cm	5.5 cm	<b>3503230</b>	351-3230SZ
250 cc	10.2 cm	4.5 cm	300 cc	10.0 cm	5.6 cm	<b>3503250</b>	
270 cc	10.4 cm	4.6 cm	325 cc	10.2 cm	5.8 cm	<b>3503270</b>	351-3270SZ
290 cc	10.8 cm	4.7 cm	350 cc	10.5 cm	5.9 cm	<b>3503290</b>	351-3290SZ
310 cc	11.0 cm	4.8 cm	375 cc	10.7 cm	6.0 cm	<b>3503310</b>	351-3310SZ
330 cc	11.3 cm	4.8 cm	400 cc	11.0 cm	6.2 cm	<b>3503330</b>	351-3330SZ
380 cc	11.7 cm	5.2 cm	450 cc	11.4 cm	6.4 cm	<b>3503380</b>	351-3380SZ
420 cc	12.0 cm	5.4 cm	500 cc	11.7 cm	6.7 cm	<b>3503420</b>	351-3420SZ
460 cc	12.4 cm	5.5 cm	550 cc	12.1 cm	6.9 cm	<b>3503460</b>	
500 cc	12.8 cm	5.6 cm	600 cc	12.4 cm	7.1 cm	<b>3503500</b>	351-3500SZ
560 cc	13.4 cm	5.7 cm	675 cc	13.1 cm	7.3 cm	<b>3503560</b>	
630 cc	13.8 cm	5.9 cm	750 cc	13.4 cm	7.4 cm	<b>3503630</b>	351-3630SZ

†Single Use Only

Individual implant dimensions may vary slightly in products of this type.  
Not all units will conform exactly to the dimensions noted above.

# SPECTRUM® POST-OPERATIVELY ADJUSTABLE SALINE BREAST IMPLANTS

## Smooth Round SPECTRUM®




Volume	NOMINAL FILL		MAXIMUM FILL	Catalog #
	Diameter	Projection	Volume	
125 cc	9.4 cm	2.7 cm	150 cc	<b>3501410</b>
175 cc	10.3 cm	3.1 cm	210 cc	<b>3501420</b>
225 cc	10.8 cm	3.5 cm	270 cc	<b>3501430</b>
275 cc	11.6 cm	3.8 cm	330 cc	<b>3501440</b>
325 cc	12.4 cm	4.0 cm	390 cc	<b>3501450</b>
375 cc	13.5 cm	4.1 cm	450 cc	<b>3501460</b>
425 cc	14.0 cm	4.2 cm	510 cc	<b>3501470</b>
475 cc	14.3 cm	4.4 cm	570 cc	<b>3501480</b>
525 cc	15.0 cm	4.4 cm	630 cc	<b>3501485</b>
575 cc	14.7 cm	4.8 cm	690 cc	<b>3501490</b>

Individual implant dimensions may vary slightly in products of this type.  
Not all units will conform exactly to the dimensions noted above.

## BREAST TISSUE EXPANDERS

### ARTOURA® PLUS Smooth and SILTEX® Breast Tissue Expanders High Profile and Ultra High Profile

All expanders include a CENTERSCOPE® Magnetic Detection Device and 21-Gauge Winged Infusion Set



High Profile				PLUS Smooth Surface 6 Suture Tabs	SILTEX® Texture Surface 3 Suture Tabs
Volume	Width	Height	Projection	Catalog #	Catalog #
225 cc	10.0 cm	10.0 cm	5.2 cm	SDC100H	TEXP100RH
300 cc	11.0 cm	11.0 cm	6.0 cm	SDC110H	TEXP110RH
375 cc	12.0 cm	12.0 cm	6.4 cm	SDC120H	TEXP120RH
475 cc	13.0 cm	13.0 cm	6.8 cm	SDC130H	TEXP130RH
500 cc	13.5 cm	13.5 cm	6.6 cm	SDC135H	TEXP135RH
600 cc	14.0 cm	14.0 cm	7.1 cm	SDC140H	TEXP140RH
750 cc	15.0 cm	15.0 cm	7.7 cm	SDC150H	TEXP150RH
850 cc	15.5 cm	15.5 cm	7.8 cm	SDC155H	TEXP155RH



Ultra High Profile				PLUS Smooth Surface 6 Suture Tabs	SILTEX® Texture Surface 3 Suture Tabs
Volume	Width	Height	Projection	Catalog #	Catalog #
350 cc	10.0 cm	10.0 cm	7.0 cm	SDC100UH	TEXP100RUH
455 cc	11.0 cm	11.0 cm	7.6 cm	SDC110UH	TEXP110RUH
535 cc	12.0 cm	12.0 cm	8.0 cm	SDC120UH	TEXP120RUH
650 cc	13.0 cm	13.0 cm	8.2 cm	SDC130UH	TEXP130RUH
700 cc	13.5 cm	13.5 cm	8.3 cm	SDC135UH	TEXP135RUH
850 cc	14.0 cm	14.0 cm	8.9 cm	SDC140UH	TEXP140RUH


PLEASE NOTE: CPX® and ARTOURA™ Breast Tissue Expanders are not MRI compatible.

Individual implant dimensions may vary slightly in products of this type.  
Not all units will conform exactly to the dimensions noted above.


## BREAST TISSUE EXPANDERS

### CPX® 4 PLUS Smooth Breast Tissue Expander with Suture Tabs


All expanders include a CENTERSCOPE® Magnetic Detection Device and 21-Gauge Winged Infusion Set



Low Height				PLUS Smooth Surface 6 Suture Tabs
Volume	Width (cm)	Height (cm)	Projection	Catalog #
250 cc	11.4 cm	8.1 cm	6.1 cm	SCPX114LH
350 cc	12.7 cm	9.4 cm	6.5 cm	SCPX127LH
450 cc	14.0 cm	10.2 cm	7.1 cm	SCPX140LH
550 cc	15.0 cm	10.9 cm	7.4 cm	SCPX150LH
650 cc	15.7 cm	11.2 cm	7.9 cm	SCPX157LH
750 cc	16.5 cm	11.9 cm	8.1 cm	SCPX165LH



Medium Height				PLUS Smooth Surface 6 Suture Tabs
Volume	Width (cm)	Height (cm)	Projection	Catalog #
275 cc	10.7 cm	9.3 cm	6.2 cm	SCPX107MH
350 cc	11.7 cm	10.0 cm	6.6 cm	SCPX117MH
450 cc	12.7 cm	10.8 cm	7.0 cm	SCPX127MH
550 cc	13.5 cm	11.7 cm	7.4 cm	SCPX135MH
650 cc	14.6 cm	12.6 cm	7.6 cm	SCPX146MH
800 cc	15.6 cm	13.3 cm	8.0 cm	SCPX156MH



Tall Height				PLUS Smooth Surface 6 Suture Tabs
Volume	Width (cm)	Height (cm)	Projection	Catalog #
250 cc	10.1 cm	10.7 cm	5.6 cm	SCPX101TH
350 cc	11.3 cm	11.8 cm	6.0 cm	SCPX113TH
450 cc	12.3 cm	12.9 cm	6.5 cm	SCPX123TH
550 cc	13.2 cm	13.8 cm	6.9 cm	SCPX132TH
650 cc	14.0 cm	14.6 cm	7.3 cm	SCPX140TH
750 cc	14.6 cm	15.3 cm	7.6 cm	SCPX146TH
850 cc	15.4 cm	15.9 cm	7.9 cm	SCPX154TH

PLEASE NOTE: CPX® Breast Tissue Expanders are not MRI compatible.

Individual implant dimensions may vary slightly in products of this type.  
Not all units will conform exactly to the dimensions noted above.

## BREAST TISSUE EXPANDERS

### CPX® 4 SILTEX® Breast Tissue Expander &

### CPX® 4 SILTEX® Breast Tissue Expander with Suture Tabs

All expanders include a CENTERSCOPE® Magnetic Detection Device and 21-Gauge Winged Infusion Set



Low Height				SILTEX® Texture Surface 3 Suture Tabs	SILTEX® Texture Surface No Suture Tabs
Volume	Width	Height	Projection	Catalog #	Catalog #
250 cc	11.4 cm	8.1 cm	6.1 cm	3549111	3548111
350 cc	12.7 cm	9.4 cm	6.5 cm	3549112	3548112
450 cc	14.0 cm	10.2 cm	7.1 cm	3549113	3548113
550 cc	15.0 cm	10.9 cm	7.4 cm	3549114	3548114
650 cc	15.7 cm	11.2 cm	7.9 cm	3549115	3548115
750 cc	16.5 cm	11.9 cm	8.1 cm	3549116	3548116



Medium Height				SILTEX® Texture Surface 3 Suture Tabs	SILTEX® Texture Surface No Suture Tabs
Volume	Width	Height	Projection	Catalog #	Catalog #
275 cc	10.7 cm	9.3 cm	6.2 cm	3549211	3548211
350 cc	11.7 cm	10.0 cm	6.6 cm	3549212	3548212
450 cc	12.7 cm	10.8 cm	7.0 cm	3549213	3548213
550 cc	13.5 cm	11.7 cm	7.4 cm	3549214	3548214
650 cc	14.6 cm	12.6 cm	7.6 cm	3549215	3548215
800 cc	15.6 cm	13.3 cm	8.0 cm	3549216	3548216



Tall Height				SILTEX® Texture Surface 3 Suture Tabs	SILTEX® Texture Surface No Suture Tabs
Volume	Width	Height	Projection	Catalog #	Catalog #
250 cc	10.1 cm	10.7 cm	5.6 cm	3549311	3548311
350 cc	11.3 cm	11.8 cm	6.0 cm	3549312	3548312
450 cc	12.3 cm	12.9 cm	6.5 cm	3549313	3548313
550 cc	13.2 cm	13.8 cm	6.9 cm	3549314	3548314
650 cc	14.0 cm	14.6 cm	7.3 cm	3549315	3548315
750 cc	14.6 cm	15.3 cm	7.6 cm	3549316	3548316
850 cc	15.4 cm	15.9 cm	7.9 cm	3549317	3548317
1445 cc	14.9 cm	16.0 cm	11.6 cm	3549336	


PLEASE NOTE: CPX® and ARTOURA™ Breast Tissue Expanders are not MRI compatible.

Individual implant dimensions may vary slightly in products of this type.  
Not all units will conform exactly to the dimensions noted above.

## TISSUE EXPANDERS

### Smooth Tissue Expanders

#### Smooth Round




Device Volume	Diameter	Projection	Catalog #
400 cc	11.3 cm	5.9 cm	<b>3504305M</b>
550 cc	13.4 cm	5.8 cm	<b>3504307M</b>
700 cc	14.6 cm	6.7 cm	<b>3504304M</b>
850 cc	15.0 cm	7.2 cm	<b>3504308M</b>

#### Smooth Elliptical




Device Volume	Length	Width	Nominal Projection	TEMPORARY OVEREXPANSION		Catalog #
				Proj.	Max. Vol.	
25 cc	6.0 cm	4.0 cm	2.3 cm	4.5 cm	75 cc	<b>3505301M</b>
50 cc	8.0 cm	5.0 cm	2.5 cm	5.3 cm	150 cc	<b>3505302M</b>
75 cc	10.0 cm	6.0 cm	2.1 cm	6.0 cm	225 cc	<b>3505303M</b>
150 cc	12.0 cm	7.0 cm	3.2 cm	7.2 cm	450 cc	<b>3505304M</b>
250 cc	14.0 cm	8.0 cm	4.5 cm	9.0 cm	750 cc	<b>3505305M</b>
125 cc	15.0 cm	5.0 cm	3.3 cm	6.5 cm	375 cc	<b>3505306M</b>

#### Smooth Crescent



Device Volume	Length	Width	Projection	Catalog #
50 cc	9.0 cm	5.7 cm	2.7 cm	<b>3504312M</b>
100 cc	11.0 cm	6.9 cm	3.7 cm	<b>3504313M</b>
200 cc	14.6 cm	7.5 cm	4.3 cm	<b>3504314M</b>
300 cc	16.5 cm	10.3 cm	4.7 cm	<b>3504315M</b>
400 cc	18.2 cm	11.4 cm	5.2 cm	<b>3504316M</b>
500 cc	19.6 cm	12.4 cm	5.5 cm	<b>3504317M</b>

#### Smooth Rectangle



Device Volume	Length	Width	Projection	Catalog #
50 cc	6.9 cm	3.2 cm	3.2 cm	<b>3504309M</b>
100 cc	6.6 cm	5.1 cm	5.6 cm	<b>3504301M</b>
250 cc	10.4 cm	7.3 cm	5.6 cm	<b>3504302M</b>
400 cc	10.6 cm	9.3 cm	6.6 cm	<b>3504303M</b>
550 cc	13.7 cm	7.1 cm	7.7 cm	<b>3504310M</b>
700 cc	15.0 cm	8.0 cm	7.9 cm	<b>3504311M</b>

Individual implant dimensions may vary slightly in products of this type.  
Not all units will conform exactly to the dimensions noted above

## ACCESSORIES

### MENTOR® Volume Sizing System

Carton contains four sterile kits. Each kit contains the following:

Description	Catalog #
MENTOR® Volume Sizing System - Display Case + Volume Sizers	<b>ESZKIT</b>
<b>Replacement Units</b>	
Display Case	<b>XXX</b>
50 cc Partial External "Gap Filling" Sizer*	<b>ESZ50</b>
75 cc Partial External "Gap Filling" Sizer*	<b>ESZ75</b>
125 cc Partial External "Gap Filling" Sizer*	<b>ESZ125</b>
150 cc External Implant Sizer*	<b>ESZ150</b>
250 cc External Implant Sizer*	<b>ESZ250</b>
350 cc External Implant Sizer*	<b>ESZ350</b>
450 cc External Implant Sizer*	<b>ESZ450</b>
550 cc External Implant Sizer*	<b>ESZ550</b>
650 cc External Implant Sizer*	<b>ESZ650</b>


### Aseptic Transfer System

Carton contains six sterile kits. Each kit contains the following:

Description	Catalog #
Aseptic Transfer System	<b>3508600</b>

Includes:

- 50 cc syringe (sterile, disposable)
- Piercing spike
- 48" tubing (saline container to syringe)
- Dual check valve
- 10" tubing (syringe to implant fill tube)
- Luer connector (attaches to implant fill tube)



### Injection Port Detector

Description	Catalog #
CENTERSCOPE® Magnetic Detection Device	<b>3504402</b>

### Injection Domes

Description	Catalog #
Dome Pack Set	<b>350DMPK</b>

### Miscellaneous

Description	Catalog #
Standard Fill Tube 8"	<b>600975001</b>
Long Fill Tube 24"	<b>600975002</b>
Check Valve	<b>200464001</b>
21-Gauge Winged Infusion Set	<b>7B3050</b>

\*A customer is able to purchase replacement units only following the purchase of the MENTOR® Volume Sizing System (ESZKIT).

PLEASE NOTE: Individual implant dimensions may vary slightly in products of this type. Not all units conform exactly to the dimensions noted above.

## The **MENTOR**<sup>®</sup> | Promise Protection Plan The Best and Most Comprehensive Implant Warranty<sup>1</sup>

MENTOR<sup>®</sup> Breast Implants are backed by the industry's most comprehensive warranty. MENTOR<sup>®</sup> Promise—with FREE and automatic comprehensive coverage—brings the peace of mind that our breast implants are protected for years to come.

---

## The **MENTOR**<sup>®</sup> Deflation Protection Program<sup>2</sup> Mentor Family of Tissue Expanders

MENTOR<sup>®</sup> is dedicated to advancing plastic surgery by listening to surgeons and customers and responding to the ever changing needs of healthcare providers in the US. The MENTOR<sup>®</sup> Deflation Protection Program is backed by our Replacement Policy for CPX<sup>®</sup> and ARTOURA<sup>™</sup> family of Tissue Expanders in the case of deflation of a MENTOR<sup>®</sup> Breast Tissue Expander in a patient.

For more information about the MENTOR<sup>®</sup> Promise Protection Plan and the MENTOR<sup>®</sup> Deflation Protection Program please visit [www.breastimplantsbymentor.com/breast-implants/warranties](http://www.breastimplantsbymentor.com/breast-implants/warranties)

### References:

1. Allergan Warranty [http://www.natrellesurgeon.com/Content/assets/Natrelle\\_Warranty%20Brochure\\_DFU](http://www.natrellesurgeon.com/Content/assets/Natrelle_Warranty%20Brochure_DFU). Accessed 04/12/2018. Sientra Product replacement policy and warranty . [http://www.sientra.com/Content/pdfs/Warranty\\_Terms\\_And\\_Conditions.pdf](http://www.sientra.com/Content/pdfs/Warranty_Terms_And_Conditions.pdf). Accessed 04/12/2018. Sientra CapCon Warranty. [http://www.sientra.com/Content/pdfs/CapCon\\_Care\\_Program\\_Terms\\_and\\_Conditions.pdf](http://www.sientra.com/Content/pdfs/CapCon_Care_Program_Terms_and_Conditions.pdf). Accessed 04/12/2018. Sientra mailing to US Plastic Surgeons
2. Mentor Replacement Policy for the CPX<sup>®</sup> and Artoura<sup>™</sup> Family of Tissue Expanders

For further information, talk to your Mentor sales representative or visit [mentorwwllc.com](http://mentorwwllc.com)



## IMPORTANT SAFETY INFORMATION:

MENTOR® MemoryGel® Breast Implants, MENTOR® MemoryShape® Breast Implants, and MENTOR® Saline-filled Breast Implants are indicated for breast augmentation in women (at least 22 years old for MemoryGel® Implants and MemoryShape® Implants, and 18 years old for Saline Implants) or for breast reconstruction. Breast implant surgery should not be performed in women with active infection anywhere in their body, with existing cancer or pre-cancer of their breast who have not received adequate treatment for those conditions, or who are currently pregnant or nursing.

Breast implants are not lifetime devices and breast implantation may not be a one-time surgery. The chance of developing complications increases over time. The most common complications for breast augmentation and reconstruction with MemoryGel® Implants include any reoperation, capsular contracture, and implant removal with or without replacement. The most common complications with MemoryShape® Implants for breast augmentation include reoperation for any reason, implant removal with or without replacement, and ptosis. The most common complications with MemoryShape® Implants for breast reconstruction include reoperation for any reason, implant removal with or without replacement, and capsular contracture. A lower risk of complication is rupture. The health consequences of a ruptured silicone gel breast implant have not been fully established. MRI screenings are recommended three years after initial implant surgery and then every two years after to detect silent rupture. Breast implants are also associated with the risk of breast implant-associated anaplastic large cell lymphoma (BIA-ALCL), an uncommon type of lymphoma. An individual's risk of developing BIA-ALCL with MENTOR® Breast Implants is low based on the incidence of worldwide cases. The most common complications with MENTOR® Saline-filled Implants include reoperation, implant removal, capsular contracture, breast pain, and implant deflation.

For MemoryGel® Implants, patients should receive a copy of Important Information for Augmentation Patients about MENTOR® MemoryGel® Breast Implants or Important Information for Reconstruction Patients about MENTOR® MemoryGel® Breast Implants. For MemoryShape® Implants, patients should receive a copy of Patient Educational Brochure - Breast Augmentation with MENTOR® MemoryShape® Breast Implants or Patient Educational Brochure - Breast Reconstruction with MENTOR® MemoryShape® Breast Implants, and a copy of Quick Facts about Breast Augmentation & Reconstruction with MENTOR® MemoryShape® Breast Implants. For MENTOR® Saline-filled Implants, patients should receive a copy of Saline-Filled Breast Implants: Making an Informed Decision. Your patient needs to read and understand the information regarding the risks and benefits of breast implants, with an opportunity to consult with you prior to deciding on surgery.

The ARTOURA® Breast Tissue Expander or CONTOUR PROFILE® Breast Tissue Expander can be utilized for breast reconstruction after mastectomy, correction of an underdeveloped breast, scar revision, and tissue defect procedures. The expander is intended for temporary subcutaneous or submuscular implantation and is not intended for use beyond six months. ARTOURA® Breast Tissue Expander and CONTOUR PROFILE® Tissue Expanders include magnetic injection domes, which contain a rare earth permanent magnet, and are NOT MRI compatible. The device could be moved by the MRI causing pain or displacement, potentially resulting in a revision surgery. Do not use the ARTOURA® Tissue Expander nor CONTOUR PROFILE® Tissue Expander in patients where an MRI may be needed. DO NOT use the ARTOURA® Breast Tissue Expander and CONTOUR PROFILE® Tissue Expander in patients that have a previously implanted device such as pacemakers, drug infusion devices, artificial sensing devices, etc. that could be affected by a magnetic field. Mentor has not tested the in vivo effects of radiation therapy with ARTOURA® Breast Tissue Expander and CONTOUR PROFILE® Expander devices and cannot warrant the safety of such use. The incidence of extrusion of the expander has been shown to increase when the expander has been placed in injured areas: scarred, heavily irradiated or burned tissue, crushed bone areas, where severe surgical reduction of the area has previously been performed and where steroids are used in the surgical pocket.

For detailed indications, contraindications, warnings, and precautions associated with the use of all MENTOR® Implantable Devices, which include MENTOR® Saline-filled Implants, MemoryGel® Implants, MemoryShape® Implants, ARTOURA® Expanders, and CONTOUR PROFILE® Expanders, please refer to the Instructions for Use (IFU) provided with each product or visit [www.mentorwllc.com](http://www.mentorwllc.com).

Important information: Prior to use, refer to the instructions for use supplied with this device for indications, contraindications, side effects, warnings and precautions.

Caution: US law restricts this device to sale by or on the order of a physician.







PART OF THE *Johnson & Johnson* FAMILY OF COMPANIES

Mentor Worldwide, LLC

Irvine, CA 92618, USA

Customer Service: (800) 235-5731 | Fax: (866) 225-2873

© Mentor Worldwide LLC 2022 020827-220114